




# CORiTEC 650i

# CORiTEC 650i/650i Loader

## WITH FUTURE-ORIENTED TECHNOLOGY

We developed the CORTEC 650i and CORTEC 650i Loader machine systems for the area of PREMIUM machines and PREMIUM requirements. The systems are offered with very attractive pricing and are equipped with high-quality industrial milling technologies, such as granite structure, linear drives in the linear axes, torque motors in the rotary axes, as well as digital absolute length measuring systems, and powerful main spindles. The machine concept impresses with its precise, vibration-free and dynamic motion sequences in the demanding and complex metal working processes. All other relevant materials can also be milled or ground with high quality using this machine system.

The CORTEC 650i loader includes a fully integrated automatic 16-fold blank changer. This enables you to operate the machine system at full capacity around the clock without supervision, while maintaining consistent high precision. Thus, the machine system is ideally suited for large labs and milling centers where these high-quality standards and large quantities are at the forefront.



- WET
- DRY
- 5 AX
- 6 mm  
32x
- Temperature control icon
- Touchscreen icon
- Video recording icon
- Blank changer icon
- 16x
- ION

# CORiTEC 650i

## HIGHLIGHTS

- Solid axis construction made of polished natural granite for 5-axis simultaneous machining in high dynamics for high-precision milling results
- absolute, high-resolution measuring systems in all axes (0.5 µm)
- high-frequency spindle up to 50.000 rpm and 3.2 kW power, with HSK 25 tool holder
- high precision due to integrated temperature compensation



WET AND DRY PROCESSING

3.2 KW HIGH-FREQUENCY SPINDLE WITH HSK-25 TOOL HOLDER

LEADING LINEAR AND TORQUE MOTOR TECHNOLOGY

POWERFUL PERFORMANCE  
high-resolution absolute-measuring systems in all axes

ERGONOMIC HANDLING  
simple touchscreen operating

MASSIVE AXLE STRUCTURE OF GRANITE NATURE

32-FOLD TOOLCHANGER

|                |                  |          |                    |                            |               |                       |                   |                  |
|----------------|------------------|----------|--------------------|----------------------------|---------------|-----------------------|-------------------|------------------|
| PMMA           | WAX              | Standard | PreMilled Abutment | Dentures                   | Crown, Bridge | Inlay, Onlay, Veneer  | Hybrid Abutment   | Splint           |
| Zr             | COMPOSITE        |          |                    |                            |               |                       |                   |                  |
| PEEK           | SINT             | C-Clamp  | Block              | zero-point clamping system | Full Denture  | Model                 | Drilling template | Model casting    |
| Glass-ceramics | PreFab Abutments |          |                    |                            |               |                       |                   |                  |
| CoCr           | Ti               |          |                    |                            | Bridge        | Telescopic technology | Abutment          | Hybrid machining |
|                |                  |          |                    |                            |               |                       |                   |                  |

# CORiTEC 650i Loader



## HIGHLIGHTS

- fully integrated automatic 16-fold blank changer
- production around the clock
- hybrid machining (Precise remilling of SLM work)
- ideally suited for large labs and milling centers
- absolute, high-resolution measuring systems in all axes (0.5 μm)
- high-frequency spindle up to 50.000 rpm and 3.2 kW power with HSK 25 tool holder
- high precision due to integrated temperature compensation

|                |                  |          |                    |                            |               |                       |                   |                  |
|----------------|------------------|----------|--------------------|----------------------------|---------------|-----------------------|-------------------|------------------|
| PMMA           | WAX              | Standard | PreMilled Abutment | Dentures                   | Crown, Bridge | Inlay, Onlay, Veneer  | Hybrid Abutment   | Splint           |
| Zr             | COMPOSITE        |          |                    |                            |               |                       |                   |                  |
| PEEK           | SINT             |          |                    |                            | Full Denture  | Model                 | Drilling template | Model casting    |
| Glass-ceramics | PreFab Abutments | C-Clamp  | Block              | zero-point clamping system |               |                       |                   |                  |
| CoCr           | Ti               |          |                    |                            | Bridge        | Telescopic technology | Abutment          | Hybrid machining |
|                |                  |          |                    |                            |               |                       |                   |                  |

# Applications & holder

## C-CLAMP

- enables 5-axis machining up to 90 degrees axis adjustment
- extension of the technical possibilities of your milling machine
- end face machining





## MANUFACTURE OF DENTAL PROSTHESES

- Baltic Denture from Merz Dental
- Dentsply Sirona Dentures



## PREFABRICATED ABUTMENTS

- short production times
- production of abutments in very simple processes
- high accuracy of fit due to prefabricated implant interface



## THERAPEUTIC SPLINTS

- bite splints, bleaching splints
- snoring splints, aligners
- simple and high-quality implementation in CAD/CAM processes




## HYBRID MACHINING - SINTERING AND MILLING

- milling post-processing of finished restorations using the SLM method
- high precision milling combined with low-cost SLM process
- compatible with most SLM systems



## GRINDING OF CAD/CAM BLOCKS

- new 6-fold adapter enables effective production
- exchangeable strips for different adapter systems
- crowns, onlays, inlays, veneers, bridge frameworks, etc.



## IMPLANT-SUPPORTED BRIDGES AND BARS / ONE-PIECE ABUTMENTS

- highest fitting accuracy due to new CAM technologies (ReFit)
- stress-free fit even with larger spans
- easy production due to high degree of automation with coordinated CAM software




## BIOHPP ELEGANCE PREFABS

- individual hybrid abutment on PEEK basis
- BioHPP for permanent dentures, free of metal, oxide and monomer
- homogeneous combination of titanium and BioHPP



## MILLING OF MODEL CASTS

- production in combustible materials, PEEK or directly in CoCr
- allergy-free and lightweight model casting when using PEEK
- cost-efficient



## MILLED IMPLANT MODELS

- DIM (Digital Implant Model) from nt-trading
- precisely positionable
- 2-part implant analog



## MODEL MILLING BAUMANN SYSTEM

- model fabrication in CAD/CAM process
- high reproducibility and precision
- model fabrication with digital impression by intraoral scanner



**T1/T6 - 3.0 | 6.0 mm**  
radius milling tool short (l = 15 mm)

|    |      |
|----|------|
| T1 | Ti   |
| T6 | CoCr |

3,0 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526029 3006

**T11 - 2.5 | 6.0 mm**  
radius milling tool (single blade, slide coated)

|     |      |
|-----|------|
| T11 | PMMA |
|     | WAX  |
|     | PEEK |

2,5 mm Ø      6 mm

up to 25 mm

Blade

Art.-No. 526004 2506

**T14/T41/T51 - 1.0 | 6.0 mm**  
radius milling tool (diamond coated)

|     |      |
|-----|------|
| T14 | Zr   |
| T41 | SINT |
| T51 | COMP |

1,0 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526013 1006

**T19 - 0.5 | 6.0 mm**  
shaft milling tool (l = 4 mm)

|     |           |
|-----|-----------|
| T19 | UNIVERSAL |
|-----|-----------|

0,5 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526001 0506

**T2/T7 - 2.0 | 6.0 mm**  
radius milling tool short (l = 12 mm)

|    |      |
|----|------|
| T2 | Ti   |
| T7 | CoCr |

2,0 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526029 2006

**T11/T13 - 2.5 | 6.0 mm**  
radius milling tool

|     |      |
|-----|------|
| T11 | PMMA |
|     | WAX  |
| T13 | Zr   |

2,5 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526019 2506

**T15/T42/T52 - 0.6 | 6.0 mm**  
radius milling tool (conical)

|     |      |
|-----|------|
| T15 | PMMA |
|     | WAX  |
|     | Zr   |
| T42 | SINT |
| T52 | COMP |

0,6 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526019 0606

**T20 - 0.6 | 6.0 mm**  
radius milling tool (conical)

|     |      |
|-----|------|
| T20 | Ti   |
|     | CoCr |

0,6 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526003 0606

**T3/T8 - 1.5 | 6.0 mm**  
radius milling tool short (l = 12 mm)

|    |      |
|----|------|
| T3 | Ti   |
| T8 | CoCr |

1,5 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526029 1506

**T12 - 1.0 | 6.0 mm**  
radius milling tool (single blade, slide coated)

|     |      |
|-----|------|
| T12 | PMMA |
|     | WAX  |
|     | PEEK |

1,0 mm Ø      6 mm

up to 25 mm

Blade

Art.-No. 526004 1006

**T16 - 1.5 | 6.0 mm**  
radius milling tool long (l = 15 mm)

|     |      |
|-----|------|
| T16 | PMMA |
|     | WAX  |
|     | PEEK |

1,5 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526019 1506

**T21 - 2.5 | 6.0 mm**  
radius milling tool (diamond)

|     |                |
|-----|----------------|
| T21 | glass ceramics |
|-----|----------------|

2,5 mm Ø      6 mm

Art.-No. 526005 2506

**T4/T9 - 1.0 | 6.0 mm**  
radius milling tool short (l = 9 mm)

|    |      |
|----|------|
| T4 | Ti   |
| T9 | CoCr |

1,0 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526029 1006

**T12/T14 - 1.0 | 6.0 mm**  
radius milling tool

|     |      |
|-----|------|
| T12 | PMMA |
|     | WAX  |
| T14 | Zr   |

1,0 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526019 1006

**T17 - 1.5 | 6.0 mm**  
shaft milling tool (l = 15 mm)

|     |      |
|-----|------|
| T17 | PMMA |
|     | WAX  |
|     | Zr   |

1,5 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526001 1506

**T22 - 1.0 | 6.0 mm**  
radius milling tool (diamond)

|     |                |
|-----|----------------|
| T22 | glass ceramics |
|-----|----------------|

1,0 mm Ø      6 mm

Art.-No. 526005 1006

**T5/T10 - 1.5 | 6.0 mm**  
shaft milling tool short (l = 12 mm, four blades)

|     |      |
|-----|------|
| T5  | Ti   |
| T10 | CoCr |

1,5 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526002 1506

**T13/T40/T50 - 2.5 | 6.0 mm**  
radius milling tool (diamond coated)

|     |      |
|-----|------|
| T13 | Zr   |
| T40 | SINT |
| T50 | COMP |

2,5 mm Ø      6 mm

up to 25 mm

Blades

Art.-No. 526013 2506

**T18 - 0.5 | 6.0 mm**  
radius milling tool (l = 4 mm)

|     |           |
|-----|-----------|
| T18 | UNIVERSAL |
|-----|-----------|

0,5 mm Ø      6 mm

up to 15 mm

Blades

Art.-No. 526026 0506

**T23 - 0.6 | 6.0 mm**  
radius milling tool (conical, diamond)


|     |                |
|-----|----------------|
| T23 | glass ceramics |
|-----|----------------|

0,6 mm Ø      6 mm

Art.-No. 526005 0606

**T26 - 3.0 | 6.0 mm**  
radius milling tool long (l = 20 mm)

|      |
|------|
| Ti   |
| CoCr |




3,0 mm Ø 6 mm

16 - 20 mm  
Blades

Art.-No. 526030 3006

**T32 - 0.6 | 6.0 mm**  
radius milling tool long (l = 12 mm)

|      |
|------|
| PMMA |
| WAX  |
| Zr   |




0,6 mm Ø 6 mm

up to 25 mm  
Blades

Art.-No. 526012 0606

**T62 - 2.0 | 6.0 mm**  
torus milling tool (r = 0,2 mm, l = 16 mm, four blades)

|      |
|------|
| Ti   |
| CoCr |




2,0 mm Ø 6 mm

up to 15 mm  
Blades

Art.-No. 526025 2006

**T80 - 2.5 | 6.0 mm**  
drilling tool (l = 22 mm)

|      |
|------|
| Ti   |
| CoCr |




2,5 mm Ø 6 mm

up to 20 mm  
Drill

Art.-No. 526023 2506

**T27 - 2.0 | 6.0 mm**  
radius milling tool long (l = 16 mm)

|      |
|------|
| Ti   |
| CoCr |




2,0 mm Ø 6 mm

16 - 20 mm  
Blades

Art.-No. 526030 2006

**T33/ T43/ T53 - 0.3 | 6.0 mm**  
radius milling tool (conical)

|      |
|------|
| PMMA |
| WAX  |
| Zr   |
| SINT |
| COMP |




0,3 mm Ø 6 mm

up to 25 mm  
Blades

Art.-No. 526019 0306

**T63 - 1.5 | 6.0 mm**  
torus milling tool (r = 0,08 mm, l = 7 mm, two blades)

|      |
|------|
| Ti   |
| CoCr |



1,5 mm Ø 6 mm

up to 15 mm  
Blades

Art.-No. 526025 1506

**T81 - 1.5 | 6.0 mm**  
drilling tool (l = 18 mm)

|      |
|------|
| Ti   |
| CoCr |




1,5 mm Ø 6 mm

up to 20 mm  
Drill

Art.-No. 526023 1506

**T28 - 1.5 | 6.0 mm**  
radius milling tool long (l = 15 mm)

|      |
|------|
| Ti   |
| CoCr |




1,5 mm Ø 6 mm

up to 25 mm  
Blades

Art.-No. 526030 1506

**T34 - 3.0 | 6.0 mm**  
radius milling tool long (l = 27 mm)

|      |
|------|
| PMMA |
| WAX  |




3,0 mm Ø 6 mm

up to 30 mm  
Blades

Art.-No. 526012 3006

**T64 - 1.5 | 6.0 mm**  
torus milling tool (r = 0,08 mm, l = 15 mm, two blades)

|      |
|------|
| Ti   |
| CoCr |




1,5 mm Ø 6 mm

up to 15 mm  
Blades

Art.-No. 526028 1506

**T98 - 2.5 | 6.0 mm**  
shaft milling tool (l = 20 mm, three blades)

|         |
|---------|
| CALIBR. |
|---------|




2,5 mm Ø 6 mm

Blades

Art.-No. 526024 2506

**T29 - 1.0 | 6.0 mm**  
radius milling tool long (l = 11 mm)

|      |
|------|
| Ti   |
| CoCr |




1,0 mm Ø 6 mm

16 - 20 mm  
Blades

Art.-No. 526030 1006

**T35 - 2.0 | 6.0 mm**  
radius milling tool long (l = 20 mm)

|      |
|------|
| PMMA |
| WAX  |




2,0 mm Ø 6 mm

up to 30 mm  
Blades

Art.-No. 526012 2006

**T67 - 2.0 | 6.0 mm**  
shaft milling tool short (l = 7 mm, four blades)

|      |
|------|
| Ti   |
| CoCr |




2,0 mm Ø 6 mm

up to 15 mm  
Blades

Art.-No. 526002 2006

**T30 - 4.0 | 6.0 mm**  
shaft milling tool long (l = 32 mm, single blade, slide coated)

|      |
|------|
| PMMA |
| WAX  |
| PEEK |




4,0 mm Ø 6 mm

up to 30 mm  
Blade

Art.-No. 526012 4006

**T61 - 3.0 | 6.0 mm**  
Quattro speed milling tool (l = 15 mm)

|      |
|------|
| CoCr |
|------|




3,0 mm Ø 6 mm

up to 15 mm  
Blades

Art.-No. 526025 3006

**T68 - 1.5 | 6.0 mm**  
shaft milling tool long (l = 15 mm, four blades)

|           |
|-----------|
| UNIVERSAL |
|-----------|



1,5 mm Ø 6 mm

16 - 20 mm  
Blades

Art.-No. 526021 1506

# CORiTEC Materials

## CORiTEC ZR (ZIRCONIUM DIOXIDE)

Covers all indications for dental restorations of up to 16 units. Ideally suited for monolithic restorations.

- highly translucent zirconium dioxide for natural esthetics
- increased strength compared to lithium dioxide

Zr ht+ Disc  
(highly translucent)




Zr Disc



Zr transpa Disc




Zr transpa Disc  
(pre-colored translucent)




## CORiTEC ZR (ZIRCONIUM DIOXIDE)

- single crowns
- single cap
- bridges up to 16 units, full bridges, and PFM technology
- good mechanical properties
- good milling properties
- corrosion resistant

CoCr Disc  
(non-precious alloy)



CoCr Mo Disc  
(non-precious alloy)



Ti Disc  
(Grad 2, Grad 4, Grad 5)




## CORiTEC MODEL DISC IVORY

- blank for the production of dental models
- low-cost material
- also suitable for test purposes or for milling calibration bodies



## JUVORA PEEK (POLYETHERETHERKETONE)

- clamp-retained constructions
- Juvora PEEK in its basic colour (brown-beige)
- highest biological compatibility



## CORiTEC WAX DISC (WAX, GRAY)

CORiTEC Wax Disc is not suitable for intraoral application. After milling, the material can be used as lost wax. Moulds can be used in casting technology. It can be milled easily and combusted residue-free.



## CORiTEC PMMA (POLYMETHYL METHACRYLATE)

- bite guards, therapeutic splints, drilling templates, bite regulators
- very good mechanical properties
- aesthetic shade effect due to high translucency
- good milling behaviour

PMMA Splint Disc



PMMA Disc







MAX PERFORMANCE  
for your milling centre

# Scanner COMPATIBLE WITH **exocad**

## CORiTEC i3Dscan eco

- all-round carefree package for newcomers or small laboratories
- simple, fast, low-priced
- very compact design



## CORiTEC i3Dscan color

- the object to be scanned is automatically guided into the measuring field
- modern touch control design
- the 180-degree opening offers a large working area



## CORiTEC i3Dscan color HR

- enormously high scanning efficiency and precision (4 µm)
- Color Texture scan
- Blue-Light LED & High resolution camera (2.8 MP)
- all modules included in the scope of delivery

# CAM software

## CORiTEC iCAM V5

- all strategies optimized for 5-axis simultaneous machining
- automatic exchange of implant connections (ReFit option)
- no annual license fees
- fully automatic and easy operation with "Wizard Workflow"



# Suction systems

## iVAC PRO+

The industrial extraction system with brushless motor, volume flow of 280 m<sup>3</sup>/hr, adjustable capacity, filter system with Teflon filter cartridge and large dust drawer, automatic compressed air self-cleaning is suitable for all dental milling machines from imes-icore®.



| iVAC PRO+              |   |
|------------------------|---|
| volume flow            | 280 m <sup>3</sup> /hr                              |
| output                 | 1200 W  |
| filter system          | Teflon filter cartridge, automatic cleaning feature |
| filter volume          | 15 liters/56.78 gpm (dust drawer)                   |
| width x depth x height | 13.8 x 13.8 x 39.4 inch / 350 x 350 x 1000 mm       |
| supply voltage         | 115 V/230 V   |
| special features       | compressed air self-cleaning system                 |

## iVAC expert

The iVAC expert is the appropriate vacuum system for the CORTEC 650i and 650i Loader milling machines. High suction power with 1000m<sup>3</sup>/hr and 100 liters filter volume and HEPA filter. The brushless turbine guarantees a long service life.




| iVAC expert            |   |
|------------------------|---|
| volume flow            | 1000 m <sup>3</sup> /hr                     |
| output                 | 2200 W                                      |
| filter system          | dust drawer, HEPA filter                    |
| filter volume          | 100 liters/378.50 gpm (filter bag)          |
| width x depth x height | 23.6 x 30.3 x 48 inch / 600 x 770 x 1220 mm |
| supply voltage         | 230 V                                       |
| special features       | high volume flow for room extraction        |

# Sintering Furnaces


## iSINT eco

The iSINT eco stands for economy and offers high-quality technology, matching accessories for high standards at a fair price. Compact and with a small footprint, the iSINT eco sintering oven still has sufficient capacity for a sintering bowl Ø 100 mm for approx. 25 units. The door stop can be mounted right and left.



## iSINT PRO

A larger sintering capacity of up to 80 single crowns is offered by the new iSINT PRO. With a heating system of four high-quality molybdenum disilicide (MoSi<sub>2</sub>) heating elements, you can choose between conventional long-term sintering or SPEED sintering at a rate of up to 99 °C/min/ 210.20°F/min. The iSINT PRO is operated with a simple and convenient program control. The programs are indicated on a 4-line LCD display. A timer function for overnight sintering or the use of pre-drying programs for shaded zirconium restorations offer further options. In addition, three service programs are available for easy maintenance of the sinter furnace.




|                             | iSINT eco           | iSINT PRO                |
|-----------------------------|---------------------|--------------------------|
| Number of heating elements  | 4                   | 4                        |
| Display                     | 7-Segment           | 4-line LCD               |
| Number of program memories  | 9                   | 30                       |
| Combustion chamber capacity | 1 x 100             | 2 x 120                  |
| max. heating rate           | 30°C/min / 86°F/min | 99 °C/min / 210.20°F/min |
| Power in W                  | 1300                | 3200                     |
| Lift function               | No                  | Yes                      |
| Speedsintering              | No                  | Yes                      |


# Special solutions Hybrid machining

## CORITEC 650i FOR SLM FINISH MACHINING


With the high-precision CORITEC 650i machine, technology partners have developed a special process to re-mill work that was previously produced in the SLM process with high precision. This allows e.g. build highly complex web constructions made of various materials such as titanium or CoCr in advance in the SLM process, and subsequently rework interfaces with the CORITEC 650i for the perfect fit and surface finish. This procedure is implemented in a special WorkFlow, in which a 3D probe is integrated in the CORITEC 650i. The touch probe determines the exact position and position of the SLM building board. Using 3D measurement functions of the CAM software, the fit and interface areas of the restorations on the SLM plate can then be precisely milled with the CORITEC 650i. This method can be used in principle with any SLM systems and with different CAM systems in combination with the CORITEC 650i.



Precise 3D measurement in the CORITEC 650i



CAM strategies for precision areas



Milling results with high surface quality

With the 3D probe, a controlled measurement is made via reference points in order to determine deviations, displacements and rotations. In the background, the imes-core® 3D-measurement software calculates the values and writes them back to the CAM software to implement the corrections.





1  
Compatible with among others the following SLM installations:



2  
Compatible with among others the following CAM systems:

# Dimensions

CORiTEC 650i



CORiTEC 650i Loader



# In numbers

CORTEC 650i

CORTEC 650i Loader

| mechanics / electronics           |   |   |
|-----------------------------------|---|---|
| number of axes and operation type | 5 axes, simultaneous machining  | 5 axes, simultaneous machining  |
| rotation angle                    | A-axis 360°, B-axis 130°  | A-axis 360°, B-axis 130°  |
| spindle / power                   | high-frequency spindle / 3.2 kW   | high-frequency spindle / 3.2 kW   |
| spindle cooling                   | water cooled  | water cooled  |
| maximum spindle speed             | 50.000 rpm  | 50.000 rpm  |
| tools                             | HSK 25  | HSK 25  |
| tool changer                      | 32-fold   | 32-fold   |
| tool length control               | ≤ 0.002 mm precision  | ≤ 0.002 mm precision  |
| blank changer                     | -   | up to 16 blanks   |
| integrated computer hardware      | Windows based   | Windows based   |
| monitor                           | integrated 15" touch screen   | integrated 15" touch screen   |
| illumination                      | 3 LED status colors   | 3 LED status colors   |
| software                          | Remote DENTAL 3.0   | Remote DENTAL 3.0   |
| air pressure                      | 6-9 bar / 87-130 PSI<br>constantly supply,<br>150 liters/min / 567.75 gpm/min | 6-9 bar / 87-130 PSI<br>constantly supply,<br>160 liters/min / 605.60 gpm/min |
| cooling liquid                    | 19 liter / 71.92 gpm integrated   | 19 liter / 71.92 gpm integrated   |
| connecting requirements           |   |   |
| weight                            | 625 kg / 1378.13 lbs  | 930 kg / 2050.65 lbs  |
| width x depth x height            | 30.9 x 43.3 x 76.4 inch<br>785 x 1100 x 1940 mm                               | 52.1 x 43.3 x 76.4 inch<br>1324 x 1100 x 1940 mm                              |
| line voltage / frequency          | 400 V - 3 stages / 50/60 Hz   | 400 V - 3 stages / 50/60 Hz   |

# Sales and support partner worldwide



Represented for you in over 100 countries.

imes-icore® GmbH has been a leading manufacturer and technology partner since 2003 in the field of dental CAD/CAM systems and solutions.

With its unique range of dental milling and grinding systems imes-icore® offers a perfect selection for all individual requirements of dental laboratories, milling centres and dental practices of all sizes.

Our “Open-System” product philosophy makes it possible to easily integrate our milling machines into your existing workflow and to integrate them with your open scanners and your CAD/CAM software. We are open for your material selection.

**imes-icore® GmbH**

Fon +49 (0) 66 72/898 228

Fax +49 (0) 66 72/898 222

dental@imes-icore.de