

Minimum footprint, maximum productivity.





The Asiga MAX is the world's most advanced digital dental 3D printer offering exceptional productivity in a small footprint. With 62µm HD print precision, the Asiga MAX is optimized for orthodontics, crown & bridge, surgical guides, dental models, custom trays, and partial dentures in lab and clinical environments. All Asiga printers are open to materials from any supplier for maximum flexibility and economy.

## Supported by Asiga's

- Lifetime technical support guarantee
- Wifi connectivity
- Automated support generation
- Touchscreen interface
- Proven manufacturing reliability

MAX		
Technical Data	MAX	MAX 🖤
PIXEL SIZE X,Y	62 µm	62 μm
BUILD SIZE X Y, Z*	119 x 67 x 75mm*	119 x 67 x 75mm*
LIGHT SOURCE	405nm	385nm
Z RESOLUTION	VARIABLE IN 1 μm	
MATERIAL SYSTEM	OPEN - USE ANY 3rd PARTY MATERIAL	
FILE INPUTS	STL, SLC, STM	
SOFTWARE	COMPOSER INCLUDED (LIFETIME SOFTWARE UPDATES INCLUDED)	
NETWORK COMPATIBILITY	WIFI AND ETHERNET	
INDUSTRY SECTORS	DIGITAL DENTISTRY / LAB PRODUCTION	
SYSTEM SIZE	260 x 380 x 370mm	
SYSTEM WEIGHT	17.5Kg	
PACKAGED SIZE/WEIGHT	410 x 500 x 480mm / 20Kg	
POWER	12VDC 10A	
PRICE RRP (USD)	\$13.990	

ouild envelope size may vary







